

Grossing Station PP1

Technical details

Workstation for 1 working person
Completely made of stainless steel
Downdraft ventilation, max. power 800m3/h
IP44 secured ventilation motor
Isolation class B
5 degree regulation of suction power
Water battery with flexible shower
(sensor operated sink battery optional)
Automatic proximity lightening
Magnetic tool holder
Shelf for paper towels
Retractable shelf for cassettes
Integrated 5l vol. canister
Perforated top desk with full cover
Integrated scale
Height adjustment legs
Power sockets for accessories (2x)
Optional Grinder

Dimensions: 1200x780x1753h



cat.no description

778100	PP1, sink on right side, no water mixer (incl.main lightening, magnetic tool holder, shelf for towels, shelf for cassettes, scale, canister)
778200	PP1, sink on right side, manual watter mixer with flexible shower (incl.main lightening, magnetic tool holder, shelf for towels, shelf for cassettes, scale, canister)
778210	PP1, sink on right side, no-touch water mixer battery powered (incl.main lightening, magnetic tool holder, shelf for towels, shelf for cassettes, scale, canister)
878800	Grinder