

Autopsy Table PSZ

Technical data

High quality stainless steel autopsy table

Cold/Hot water flexible hand shower

Superb stability

Sink dimensions 450x450x200mm

Optional dowdraft fume extraction

Optional height adjustment

Optional body trays


Optional head rest

Optional organ table

Optional tool table

dimensions: 2700x850x750h fixed version

dimensions: 2700x850x750-1000h elevated version



cat.no.: description

PSZ-P	Fixed autopsy table
PSZ-O	Fixed autopsy table with fume extraction system
PSZ-ZO	Elevated autopsy table with fume extraction system
PSZ-Z	Elevated autopsy table

* no organ and no tools table included