


High-capacity archive for slide storage

Z56000



- Storage cabinet for slides storage
- Complete steel construction
- Easy to manipulate
- 140 shelves for easy transport to your desk
- dimensions: 535x510x1140mm
- Capacity 56.000 slides
- Shelves are removable



Elegant blu-gray design for your lab

Great stability

Space savings

Z26000



- Slide storage archiving cabinet
- 10 independent removable drawers
- Each drawer with cover
- Complete steel construction
- Dimensions: 305x395x1170mm
- Elegant gray color
- Easy transport of drawers to your desk
- Capacity 26.000 slides



Elegant gray design for your lab

Great stability

Space savings

Z8000



- Storage cabinet for paraffin block
- 15 drawers with possibility to remove
- Blue-gray styled color
- Complete steel construction
- Dimensions: 450x350x940mm
- Possible mounting second cabinet on top of the first one
- Easy drawer transport to your desk
- Capacity 8.000 blocks.

